


## Product Datasheet

### ZNF592 antibody (orb227156)

|                            |                                      |
|----------------------------|--------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to ZNF592 |
| <b>Species/Host</b>        | Rabbit                               |
| <b>Reactivity</b>          | Human, Mouse                         |
| <b>Conjugation</b>         | Unconjugated                         |
| <b>Tested Applications</b> | IHC-P, WB                            |
| <b>Note</b>                | For research use only                |
| <b>Clonality</b>           | Polyclonal                           |
| <b>Dilution Range</b>      | WB: 1:500-1:1000, IHC-P: 1:50-1:200  |
| <b>Expiration Date</b>     | 12 months from date of receipt.      |



Western blot analysis of K562 cells using ZNF592 antibody.