


Product Datasheet

Histone H4 (Acetyl-Lys5) antibody (orb227349)

Description	Rabbit polyclonal antibody to Histone H4 (Acetyl-Lys5)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200, IF/ICC: 1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of COS7 cells using...