


## Product Datasheet

### EDG5 antibody (orb227378)

|                            |                                    |
|----------------------------|------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to EDG5 |
| <b>Species/Host</b>        | Rabbit                             |
| <b>Reactivity</b>          | Human, Mouse, Rat                  |
| <b>Conjugation</b>         | Unconjugated                       |
| <b>Tested Applications</b> | WB                                 |
| <b>Note</b>                | For research use only              |
| <b>Clonality</b>           | Polyclonal                         |
| <b>Dilution Range</b>      | WB: 1:500-1:1000                   |
| <b>Expiration Date</b>     | 12 months from date of receipt.    |



Western blot analysis of COLO205 cells u...