


Product Datasheet

GPR132 antibody (orb227391)

Description	Rabbit polyclonal antibody to GPR132
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:1000, IF/ICC: 1:100-1:500, ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of COS-7 cells using GPR132 antibody.



IF image of hela cells lysates using GPR132 antibody.