


Product Datasheet

GPR158 antibody (orb227454)

Description	Rabbit polyclonal antibody to GPR158 also known as probable G-protein coupled receptor 158 whose exp...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1000, IHC-P: 1:50-200, IF/ICC: 1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 cells using...