


## Product Datasheet

### LGR6 antibody (orb227476)

|                            |                                       |
|----------------------------|---------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to LGR6    |
| <b>Species/Host</b>        | Rabbit                                |
| <b>Reactivity</b>          | Human                                 |
| <b>Conjugation</b>         | Unconjugated                          |
| <b>Tested Applications</b> | ICC, IF, WB                           |
| <b>Note</b>                | For research use only                 |
| <b>Clonality</b>           | Polyclonal                            |
| <b>Dilution Range</b>      | WB: 1:500-1:1000, IF/ICC: 1:100-1:500 |
| <b>Expiration Date</b>     | 12 months from date of receipt.       |



Western blot analysis of HUVEC cells using LGR6 antibody (orb227476).