


Product Datasheet

TAS2R14 antibody (orb227649)

Description	Rabbit polyclonal antibody to TAS2R14
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:1000, IF/ICC: 1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF-7 cells using TAS2R14 antibody (orb227649).