


Product Datasheet

Integrin Alpha3 (CD49c) antibody (orb227721)

Description	Rabbit polyclonal antibody to Integrin Alpha3 (CD49c)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1-500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cell lysate using Integrin Alpha3 (CD49c) antibody (orb227721).