


Product Datasheet

RUNX2 antibody (orb227722)

Description	Rabbit polyclonal antibody to RUNX2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200, IF/ICC: 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using RUNX2 antibody (orb227722).