


Product Datasheet

Interleukin 5 antibody (orb227775)

Description	Rabbit polyclonal antibody to Interleukin 5
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HT-29 cells using Interleukin 5 antibody (orb227775).