


Product Datasheet

Annexin A3 antibody (orb227810)

Description	Rabbit polyclonal antibody to Annexin A3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cells using Annexin A3 antibody (orb227810).