


Product Datasheet

ABCG2 antibody (orb227828)

Description	Rabbit polyclonal antibody to ABCG2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HT-29 cells using ABCG2 antibody (orb227828).