


## Product Datasheet

### MDR1 antibody (orb227836)

|                            |                                    |
|----------------------------|------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to MDR1 |
| <b>Species/Host</b>        | Rabbit                             |
| <b>Reactivity</b>          | Human, Mouse, Rat                  |
| <b>Conjugation</b>         | Unconjugated                       |
| <b>Tested Applications</b> | IHC-P                              |
| <b>Note</b>                | For research use only              |
| <b>Clonality</b>           | Polyclonal                         |
| <b>Dilution Range</b>      | IHC-P: 1:50-1:200                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.    |



Immunohistochemical  
staining of human  
br...