


Product Datasheet

MMP-9 antibody (orb227878)

Description	Rabbit polyclonal antibody to MMP-9
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-P: 1:50-1:200-P WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis on HepG2 cell lysates using MMP-9 antibody. The lane on the left is treated with the antigen-specific peptide.

Western blot analysis of HepG2 cells using...