


Product Datasheet

Caspase 9 (Cleaved-Asp353) antibody (orb227889)

Description	Rabbit polyclonal antibody to Caspase 9 (Cleaved-Asp353)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH/3T3 cells u...