


Product Datasheet

Collagen IV Alpha2 antibody (orb227913)

Description	Rabbit polyclonal antibody to Collagen IV Alpha2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF/ICC: 1:100-1:500, ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of Hela cells lysa...



Immunohistochemical staining of human Br...



Immunofluorescent analysis of NIH-3T3 ce...