



Product Datasheet

Notch 1 (Cleaved-Val1744) antibody (orb227956)

Description	Rabbit polyclonal antibody to Notch 1 (Cleaved-Val1744)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IHC-P 1:50-1:200, ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse sp...



Western blot analysis of Mouse brain usi...