


## Product Datasheet

### Caveolin 2 (Phospho-Tyr27) antibody (orb228068)

|                            |                                          |
|----------------------------|------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Caveolin 2 |
| <b>Species/Host</b>        | Rabbit                                   |
| <b>Reactivity</b>          | Human                                    |
| <b>Conjugation</b>         | Unconjugated                             |
| <b>Tested Applications</b> | WB                                       |
| <b>Note</b>                | For research use only                    |
| <b>Clonality</b>           | Polyclonal                               |
| <b>Dilution Range</b>      | WB: 1:500-1:2000                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.          |



Western blot analysis of human fetal lung tissue.