


Product Datasheet

Collagen III antibody (orb228096)

Description	Rabbit polyclonal antibody to Collagen III
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-P: 1:50-1:200, IF/ICC: 1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HT1080 Whole ce...