


Product Datasheet

Alpha 1 Fetoprotein antibody (orb228142)

Description	Rabbit polyclonal antibody to Alpha 1 Fetoprotein
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF/ICC: 1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of various cell li...