



Product Datasheet

VCAM1 antibody (orb228217)

Description	Rabbit polyclonal antibody to VCAM1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500 - 2000, IHC-P: 1:50 - 200, ELISA: 1:20000 - 40000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human
Ly...



Western blot analysis
of extracts from r...