


Product Datasheet

VHL antibody (orb228239)

Description	Rabbit polyclonal antibody to VHL also known as protein G7 whose expression can be detected in adult...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of various cell li...