



Product Datasheet

ACVR1C antibody (orb228253)

Description	Rabbit polyclonal antibody to ACVR1C
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, P-ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis
of 293T cell lysat...



Immunohistochemical
staining of human
br...