



Product Datasheet

NeuN antibody (orb228280)

Description	Rabbit polyclonal antibody to NeuN
Species/Host	Rabbit
Reactivity	Bovine, Rabbit, Sheep
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, ELISA: 1:20000-40000
Expiration Date	12 months from date of receipt.



Western blot analysis of various cell li...



Immunohistochemical staining of Mouse br...