


## Product Datasheet

### Osteopontin antibody (orb228530)

|                            |                                           |
|----------------------------|-------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Osteopontin |
| <b>Species/Host</b>        | Rabbit                                    |
| <b>Reactivity</b>          | Human, Mouse, Rat                         |
| <b>Conjugation</b>         | Unconjugated                              |
| <b>Tested Applications</b> | IHC-P, WB                                 |
| <b>Note</b>                | For research use only                     |
| <b>Clonality</b>           | Polyclonal                                |
| <b>Dilution Range</b>      | WB: 1:500-1:2000, IHC-P: 1:50-1:200       |
| <b>Expiration Date</b>     | 12 months from date of receipt.           |



Western blot analysis of HeLa cell lines...