



Product Datasheet

NR1I2 antibody (orb228613)


Description	Rabbit polyclonal antibody to PXR
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF: 1:100-500, P-ELISA: 1:20000-40000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse liver usi...



Immunofluorescence analysis of HeLa cell...



Immunohistochemical staining of paraffin...