



Product Datasheet

PDL1 antibody (orb228661)

Description	Rabbit polyclonal antibody to PDL1
Species/Host	Rabbit
Reactivity	Bovine, Human, Mouse, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-2000, ELISA: 1:20000-40000
Expiration Date	12 months from date of receipt.



Western blot analysis
of rat lung, MCF7 ...



Immunohistochemical
staining of human
CO...