



Product Datasheet

SAA1 antibody (orb228668)

Description	Rabbit polyclonal antibody to SAA1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, ELISA(peptide): 1:20000-1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of human liver cel...



Immunohistochemical staining of Mouse li...