


Product Datasheet

GRIA1 antibody (orb228776)

Description	Rabbit polyclonal antibody to GRIA1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:500-1:2000, IF/ICC: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of rat brain cell

...