



Product Datasheet

KLHL3 antibody (orb2288)


Description	KLHL3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human KLHL3 (51-160/587 aa)
Target	KLHL3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65 kDa
Uniprot ID	Q9UHH7
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of human HL-60 lys...



Western blot analysis of human A549 lysa...



Immunohistochemical staining of Rat kidn...