



Product Datasheet

CD70 antibody (orb228901)


Description	Rabbit polyclonal antibody to CD70
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:100-1:500, ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of various cell li...



IF analysis of HuVEC cells lysate using ...



IHC-P image of human pancreas tissue usi...