


Product Datasheet

DNMT3L antibody (orb229063)

Description	Rabbit polyclonal antibody to DNMT3L
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ChIP, ICC, IF, IHC-P, IP, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF/ICC: 1:20-1:50, IP: 1:20-1:50, ChIP: 1:20-1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of various cell li...