


Product Datasheet

Histone H3K4me3 antibody (orb229070)

Description	Rabbit polyclonal antibody to Histone H3K4me3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, ICC, IF, IHC-P, IP, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF/ICC: 1:50-1:200, IP: 1:50-1:200, ChIP: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cell line ...