


Product Datasheet

H3K9me3 antibody (orb229073)

Description	Rabbit polyclonal antibody to H3K9me3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, ICC, IF, IHC-P, IP, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200, IP: 1:50-200, ChIP: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cell line

...