


Product Datasheet

Histone H3K79me2 antibody (orb229081)

Description	Rabbit polyclonal antibody to Histone H3K79me2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, ICC, IF, IHC-P, IP, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:100, IF/ICC: 1:50-1:200, IP: 1:50-1:200, ChIP: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cell line ...