



Product Datasheet

IL8 antibody (orb229133)


Description	Rabbit polyclonal antibody to interleukin 8 is also known as IL8/NAP1 form V and MDNCF whose express...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from h...



Effect of JHQ on core target proteins ab...



Effects of inhibition of HIF-1α on expe...