


Product Datasheet

DNMT3B antibody (orb229237)

Description	Rabbit polyclonal antibody to DNMT3B
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:200-1:1000, IF/ICC: 1:100-1:500, ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 cell line ...



Immunofluorescence analysis of HepG2 cel...