



Product Datasheet

GRM1 antibody (orb229255)

Description	Rabbit polyclonal antibody to GRM1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100, ELISA: 1:20000-40000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cell lysat...



Immunohistochemical staining of Human li...