


## Product Datasheet

### GRK1 antibody (orb229257)

|                            |                                     |
|----------------------------|-------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to GRK1  |
| <b>Species/Host</b>        | Rabbit                              |
| <b>Reactivity</b>          | Human, Mouse, Rat                   |
| <b>Conjugation</b>         | Unconjugated                        |
| <b>Tested Applications</b> | IHC-P, WB                           |
| <b>Note</b>                | For research use only               |
| <b>Clonality</b>           | Polyclonal                          |
| <b>Dilution Range</b>      | WB: 1:500-1:2000, IHC-P: 1:50-1:200 |
| <b>Expiration Date</b>     | 12 months from date of receipt.     |



Western blot analysis of Mouse spleentis...