


Product Datasheet

SLC16A1 antibody (orb229279)

Description	Rabbit polyclonal antibody to SLC16A1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot
analysis of
Mouseseminoma
t...