


Product Datasheet

H3K14me3 antibody (orb229380)

Description	Rabbit polyclonal antibody to H3K14me3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF/ICC: 1:20-1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa and H3 pro...