


Product Datasheet

MLKL antibody (orb229539)

Description	Rabbit polyclonal antibody to MLKL also known as hMLKL whose expression can be detected in the adult...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of THP-1 cell line...