



Product Datasheet

C5orf51 antibody (orb22954)


Description	C5orf51 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human UPF0600 protein C5orf51 protein (8-152AA)
Target	C5orf51
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A6NDU8
Dilution Range	IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human th...



Immunohistochemical staining of human te...



Immunofluorescence analysis of HepG-2 ce...