


Product Datasheet

Retinoblastoma (Phospho-Thr826) antibody (orb229603)

Description	Rabbit polyclonal antibody to Retinoblastoma
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IF/ICC: 1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of MOLT cell lysat...