



Product Datasheet

AKT (Phospho-Ser473) antibody (orb229611)

Description	Rabbit polyclonal antibody to AKT (Phospho-Ser473)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF/ICC: 1:100-1:500, ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human lu...



Immunofluorescence analysis of 293 cell ...



Immunofluorescence analysis of 293 cell ...