


## Product Datasheet

### AKT (Phospho-Ser473) antibody (orb229612)

|                            |                                                          |
|----------------------------|----------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to AKT (Phospho-Ser473)       |
| <b>Species/Host</b>        | Rabbit                                                   |
| <b>Reactivity</b>          | Human, Mouse, Rat                                        |
| <b>Conjugation</b>         | Unconjugated                                             |
| <b>Tested Applications</b> | ICC, IF, IHC-P, WB                                       |
| <b>Note</b>                | For research use only                                    |
| <b>Clonality</b>           | Polyclonal                                               |
| <b>Dilution Range</b>      | WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF/ICC: 1:100-1:500 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                          |



Western blot analysis of 293 cell lysate...