


## Product Datasheet

### 2A6 (Cytochrome P450) antibody (orb229668)

|                            |                                   |
|----------------------------|-----------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to 2A6 |
| <b>Species/Host</b>        | Rabbit                            |
| <b>Reactivity</b>          | Human                             |
| <b>Conjugation</b>         | Unconjugated                      |
| <b>Tested Applications</b> | WB                                |
| <b>Note</b>                | For research use only             |
| <b>Clonality</b>           | Polyclonal                        |
| <b>Dilution Range</b>      | WB: 1:500-1:1000                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |



Western blot analysis of Jurkat cells us...