


## Product Datasheet

### Cytochrome P450 2B1 antibody (orb229670)

|                            |                                                   |
|----------------------------|---------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Cytochrome P450 2B1 |
| <b>Species/Host</b>        | Rabbit                                            |
| <b>Reactivity</b>          | Human                                             |
| <b>Conjugation</b>         | Unconjugated                                      |
| <b>Tested Applications</b> | WB                                                |
| <b>Note</b>                | For research use only                             |
| <b>Clonality</b>           | Polyclonal                                        |
| <b>Dilution Range</b>      | WB: 1:500-1:1000                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                   |



Western blot analysis of HT-29 cells using...