



Product Datasheet

GAB1 Antibody (orb229776)


Description	Mouse monoclonal antibody to GAB1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Purified recombinant fragment of human GAB1 (AA: 661-724) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	32.4kDa
Entrez	2549
Dilution Range	RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



ELISA Black line: Control
Antigen (100 ng...)



Flow cytometric analysis of HeLa cells u...



Western blot analysis of human recombina...