



Product Datasheet

LDLR Antibody (orb229780)


Description	Mouse monoclonal antibody to LDLR
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC
Immunogen	Purified recombinant fragment of human LDLR (AA: 22-150) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	95.4kDa
Entrez	3949
Dilution Range	IHC: 1/200 - 1/1000; RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



ELISA Black line:
Control Antigen (100
n...



Flow cytometric
analysis of HeLa cells
U...



Immunohistochemical
analysis of formalin...